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 serf

 serf
EXTREMITIES




<https://doc.serf.fr/0916.pdf>

- FR** Notice d'instructions
Instruments réutilisables SERF ... **13**
- EN** Instructions for use
SERF reusable instruments ... **17**
- AR** دليل التعليمات
سيرف أدوات يمكن إعادة استعماله ... **21**
- BG** Инструкция за употреба
Нструменти за многократна употреба SERF ... **24**
- DA** Brugsanvisning
Genanvendelige instrumenter fra SERF ... **29**
- DE** Gebrauchsanweisung
Wiederverwendbare Instrumente von SERF ... **33**
- ES** Instrucciones de uso
Instrumentos reutilizables SERF ... **38**
- FI** Käyttöohje
Uudelleen käytettävät instrumentit SERF ... **42**
- HU** Tájékoztató
SERF többször használatos műszerek ... **46**
- IT** Avvertenza di istruzioni
Strumenti riutilizzabili SERF ... **50**
- KA** ინსტრუქციის ფურცელი
მრავალჯერადი ხელსაწყოები SERF ... **54**
- LT** Naudojimo instrukcija
SERF daugkartinio naudojimo instrumentai ... **59**
- NB** Bruksanvisning
Gjenbrukbare instrumenter fra SERF ... **63**
- NL** Gebruiksaanwijzing
Herbruikbare instrumenten SERF ... **67**
- PL** Instrukcja obsługi
Instrumenty wielokrotnego użytku SERF ... **71**
- PT** Instruções de utilização
Instrumentos reutilizáveis SERF ... **75**
- RU** Инструкция
Инструменты многоразового использования SERF ... **79**
- SK** Návod na použitie
Nástroje opakovane použiteľné SERF ... **84**
- SV** Instruktionsblad
Återanvändbara SERF-instrument ... **88**
- VI** Hướng dẫn sử dụng
Công cụ tái sử dụng SERF ... **92**

1

2

3

4


5




6

7

8



9




10

11

12



YES**NO**

ICJB E

REF RA90430072

14

Components:

BAJB2

REF RA90430161

COJB 12




REF RA90430011


EICJB


REF RA90250010

MV 700

REF RA90450030


MP020L

REF RA90460171

1. VERSIONS

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1. VERSIONS

Date	Status / Nature of modifications	URL Internet
June 2024	CE (MDR) certification of "hip" baskets (section "14. Date instruments first put on the market").	https://doc.serf.fr/0916.pdf
April 2024	- Integration of instruments assembly and disassembly instructions (updated instructions for use content accordingly). - Addition of Lithuanian and Polish languages. - Deletion of the Chinese version of the texts in the notice.	https://doc.serf.fr/0916-rev5-1.pdf
October 2022	- New corporate charter (logotype, model). - Addition of a new section "3. Extremity instrument distributor". - Addition of a "Distributor" symbol (section "5. Packaging"). - Addition of a clarification on the life of instruments (section "7. Reprocessing limits of instruments"). - Modification of the paragraph "Pre-cleaning preparation" (section "8. Reprocessing instructions"). - Modification of the section "14. Date instruments first put on the market".	https://doc.serf.fr/0916-ed5.pdf
August 2021	- Addition "SERF Extremity" logotype (cover). - Distinction of devices in compliance with regulation 2017/745 and devices according to directive 93/42/EEC. - Specification of scope of the instrumentation supported by these instructions (reusable and non-sterile instruments). - Addition symbol "Date of manufacture" (section "Packaging"). - Addition of an item concerning the notification of serious accidents (section "Warnings"). - Modification sterilization parameters (section "Reprocessing instructions"). - Addition of Arabic language.	https://doc.serf.fr/0916-ed4.pdf
January 2021	Overall modification of the instruction leaflet.	https://doc.serf.fr/0916-ed3.pdf
December 2019	Modification section "Reprocessing instructions".	https://doc.serf.fr/0916-ed2.pdf
September 2019	Creation	https://doc.serf.fr/0916-ed1.pdf

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3. EXTREMITY INSTRUMENT DISTRIBUTOR

SERF Extremity
85 avenue des Bruyères
69150 Décines-Charpieu - France
Tel. +33 (0)4 72 05 60 26
Fax +33 (0)4 72 02 19 18
Mail: serf.extremity@serf.fr
Internet: www.serf.fr

4. GENERAL INFORMATION ON THE INSTRUMENTS

- This document applies to all the reusable instruments manufactured by SERF and delivered non-sterile.
- The non-sterile, reusable instruments are used to insert the associated SERF implants.
- Each instrument kit is described in the accompanying implantation sheet.
- For sterile or single-use instruments, please refer to the specific instructions.

5. PACKAGING

The meanings of the regulatory symbols and pictograms on the instrument packaging labels are given below:

	Catalogue reference
	Batch code
	Unique device identification
	Device ID
	Production unit
	Medical device
	Number of devices in the package
	CE marking (European Conformity)
	CE marking, followed by number identifying the notified body
	Date of manufacture
	Manufacturer
	Non sterile

See the instructions for use

Distributor

6. WARNINGS

- These instructions are provided to guide the medical staff during the reprocessing phase of surgical instruments. These instructions have been approved by SERF. It is the medical institution's responsibility to:
 - ensure that the reprocessing is performed using the required equipment,
 - and that all the steps are carried out by trained staff.
- Non-sterile and reusable instruments must be reprocessed (cleaned, inspected and sterilised) in compliance with the instructions described hereafter prior to use.
- Be careful with sharp instruments.
- Any serious instrument-related incident must be notified to SERF and the relevant authority of the member State in which the user and/or the patient is established.
- The instruments must be cleaned out of the sterilisation trays. The trays are designed for the sterilisation and storage of the instruments.

- **Composition of the detergent and/or disinfectant products:** using a hydrogen peroxide- (H_2O_2), nitric acid- (HNO_3), sulfuric- (H_2SO_4) or oxalic-based solution is not recommended for cleaning instruments made of titanium alloy, or aluminium boxes (may turn them black). Similarly, we do not authorise any immersion of our instruments in solutions containing aldehyde or a strong concentration of chlorine (see sodium hypochlorite).
- **Location of instruments or fragments of them:** instruments or fragments of instruments can be located visually without necessarily having to use an external imaging device. As metallic instruments or fragments of them are radio-opaque, they can however be located using radiography according to their size.

7. REPROCESSING LIMITS OF INSTRUMENTS

The successive cleaning and sterilization cycles (in compliance with the instructions in this instruction booklet) have undesirable effects on the lifetime and function of the instruments. Their life expectancy is estimated at 81 uses, equivalent to a 3-year useful life. The wear and damage can be materialised by corrosion, discoloration, scratches or impacts (non-exhaustive list). Instruments with defects of function or identification must not be used. A photo library is available for this purpose in this instruction booklet; it illustrates the limits of use of instruments (see section "12. Limits of use of instruments" and figures ① to ⑬).

8. REPROCESSING INSTRUCTIONS

- The following instructions allow medical devices to be prepared with a view to being reused.
- The description of the products, materials and equipment, as well as their conditions of use specified in each of the stages, will encourage the safe reprocessing of SERF instruments.
- It is the surgeon's responsibility to always ensure that the personnel, the equipment and the process involved in reprocessing SERF devices achieve the desired result. This generally requires a regular validation and inspection of the processes. Any deviation from the given instructions by the relevant department should also be correctly assessed, to determine the effectiveness or any adverse effects.
- The maximum period of contact between each instrument and the patient's body is 10 minutes.

- The products, dosages, equipment and conditions of use given in brackets are part of the SERF protocol for preparing equipment for sterilisation. The product references quoted have been tested by SERF and are given as examples.
- Prior to any use, it is essential to clean, decontaminate, inspect and sterilise the instruments, respecting the instructions below.

Places of use

Important: after use, clean and decontaminate immediately at the end of the operation.

It is essential to prevent dirt drying on the equipment and all recyclable equipment must be put to soak straight away after each use. Moreover, the largest stains (blood, tissue, debris) should be removed beforehand using lint-free sponges or brushes, or by rinsing with cold water. Prolonged contact with blood, bones fragments and other residues may harm the instruments and compromise the effectiveness of the treatment overall.

Do not put dirty instruments back in their storage box.

The instruments must be disassembled to be cleaned in compliance with the instructions provided with the instrument kit. When no disassembly instructions are available, the instrument must not be disassembled for cleaning.

Confinement and transport

After surgery, immerse the instruments in a neutral aqueous solution (e.g. physiological liquid or distilled water) for 15 minutes.

The surgical instruments used are considered to present a biological risk and must be taken to a specific reprocessing zone in compliance with the medical institution's policy.

Pre-cleaning preparation

It is recommended that instruments be reprocessed immediately after use.

The instruments must be cleaned out of the trays and lids used for their storage. Particular attention must be paid to sharp instruments. It is recommended that these be cleaned separately to reduce the risks of injury.

The instruments must be disassembled for cleaning in compliance with the instructions presented in section "13. Instructions for assembly and disassembly of instruments". When no disassembly instructions are

supplied, the instrument must not be disassembled for cleaning.

Any instruments with nozzles must be purged to remove debris and brushed using a brush that is appropriate for the nozzle size.

Automated cleaning

1. Place the instruments in the washer-disinfector, then make sure that the instruments and nozzles have maximum exposure to the detergent and rinsing water and that it can flow freely.

2. The instruments must be disassembled for cleaning in compliance with the instructions presented in section "13. Instructions for assembly and disassembly of instruments". When no disassembly instructions are available, the instrument must not be disassembled for cleaning.

3. Automated cleaning must be carried out in a certified washer-disinfector in compliance with ISO 15883-1 or an equivalent standard. The following instructions have been approved by the manufacturer:

Phase	Duration (mn)	Water temperature	Detergent / Water quality
Pre-wash	2	< 40°C	Raw water
Washing	10	55°C	Mediclean 0.5% Raw water
Rinsing	3	Ambient	Raw water
Rinsing	2	Ambient	Reverse osmosis water

Automated disinfection

Thermal disinfection is recommended to secure the handling of the instruments prior to steam sterilisation.

Thermal disinfection must be carried out in a washer-disinfector in compliance with ISO 15883-1 and ISO 15883-2, or an equivalent standard. The thermal disinfection in the washer-disinfector must be approved to obtain an A_0 of 3000.

Place the instruments in such a way as to ensure that all the water drains away, cannulated devices must be placed vertically.

Phase	Duration (mn)	Water temperature	Detergent / Water quality
Thermal disinfection	6	93°C	Mediklar 0,05% Reverse osmosis water

Drying

Drying must be done in a washer-disinfector in compliance with ISO 15883-1 and ISO 15883-2, or an equivalent standard. The effectiveness of drying in a washer-disinfector can vary considerably according to its design and load.

The following cycle is approved by the manufacturer:

Phase	Duration (mn)	Air temperature	Air type
Drying	30	110°C	Medical grade

If you are using a washer/disinfector, temperature must not exceed 120°C, (air drying 110°C, then air 30°C for cooling).

After automated drying, inspect the instruments. If there are traces of moisture on an instrument, it must be dried manually.

Manual drying should be done using a clean, absorbent nonwoven support or a clean, lint-free cloth together with "medical" compressed air.

Inspection and testing

Inspect all the instruments visually to check that there is no dirt, damage or moisture.

Inspect all the instruments to check:

- **the presence of any damage** such as corrosion, discoloration, scratches, impacts or wear,
- **reassembly:** inspect the instruments that were disassembled and check that they are assembled in accordance with the instructions presented in section "13. Instructions for assembly and disassembly of instruments",
- **absence of moisture:** carefully inspect devices with nozzles and moving components; if moisture is detected, dry the device manually in accordance with the drying instructions given earlier in the section "Drying",
- **cleanliness:** if dirt is identified, repeat the cleaning steps until all the dirt is removed from the device,
- **the instrument's functionality,** for example: the instruments' cutting capacity, bent axes, movement of hinges, joints or the locking system, moving components such as handles, lock mechanisms, etc.,

- **check that all of the instruments are present** and that they are placed in the trays as shown in the instructions sent with each instrument kit.

Instruments with worn marking, instruments that cannot be identified and damaged or non-functional instruments must be discarded.

Illustrations of the acceptance criteria of instruments during the inspection stage are given in figures ❶ to ❸ and section "12. Limits of use of instruments".

Packaging

Note: instruments must be packaged rapidly to prevent any recontamination.

Use only legally approved, medical quality sterile barriers (bag, wrap, container) for the final packaging of steam-sterilised instruments.

Individual packaging: use an appropriate material for steam sterilisation and one that complies with ISO 11607-1. Ensure that the bag is large enough to take the instrument without stretching at the edges. Do not put several instruments in the same bag.

Batch packaging: the instruments must be replaced in their original kit as shown in the instructions sent with each instrument kit. No additional instruments should be placed in the tray. Check that the instruments are not touching, particularly instruments with sharp edges. Package the tray using one of the following appropriate methods:

- in an approved sterilisation container with a watertight lid,
- in steam sterilisation packaging of medical quality standard, using the AAMI double-wrapping method or equivalent.

Sterilisation

Sterilisation in a steam pressure autoclave must be done according to a locally approved vacuum cycle. The autoclave must be validated to local requirements and to ISO 17665-1. The autoclave must be installed and maintained in compliance with local standards and the manufacturer's requirements.

Use an autoclave cycle designed to evacuate porous or cannulated devices according to the manufacturer's instructions. Do not exceed the autoclave's maximum load.

The following cycles are validated for a Sterility Assurance Level (SAL) of 10^{-6} :

Cycle type	Minimum exposure time	Temperature	Drying time
Pre-vacuum	4 mn	132°C	30 mn
Pre-vacuum	3 mn	134°C	30 mn

However, a longer exposure time may be used in compliance with local requirements, like, for example the 134°C cycle for 18 minutes.

Storage

The sterile instruments must be stored in a dry, clean environment, protected from direct light, pests and extreme temperatures and humidity.

Refer to suppliers' instructions for double wrapping or containers relative to time limits of storage and requirements of temperature and humidity for storage.

The post-sterilization instrument packaging must be carefully examined before being opened, to check packaging integrity.

9. INACTIVATION OF NONCONVENTIONAL TRANSMISSIBLE AGENTS (NCTA) IF REQUIRED

- After the cleaning phase, any equipment that requires it must undergo an NCTA inactivation process. The client must choose the inactivation process in accordance with instruction DGS/R13/2011/449 dated December 1st 2011, published by the French Ministry of Health; instruction relative to updating recommendations aimed at reducing the risks of NCTA transmission during invasive procedures. We recommend immersing our instruments in soda (optimum concentration is 1 mol/l of sodium hydroxide/NaOH) for 1 hour, then steam sterilisation at 134°C for 18 minutes minimum in a porous load autoclave (pressurised steam steriliser). We do not allow the immersion of our ancillary equipment in sodium hypochlorite. It is, moreover, formally forbidden to submerge our boxes.
- NCTA inactivation must always be followed by careful rinsing in water. Equipment must be dried before packing.

10. RETURNING LOANED EQUIPMENT TO SERF

- Following the treatments listed above, put the instruments back carefully in their respective places in the kit boxes and containers in such a way as to prevent any damage during transport, then fill

in the "traceability/treatment" form (or exchange form if deposit), then send the cleaned and decontaminated equipment back to SERF.

- **Before returning loaned equipment to SERF**, inspect the instruments that were disassembled and check that they are reassembled in accordance with the instructions presented in section "13. Instructions for assembly and disassembly of instruments".

11. DISPOSAL OF INSTRUMENTS

Sorting of waste with infectious risk and waste similar to household waste is the responsibility of the professional who has produced the waste.

Generally speaking, ordinary waste such as packaging can be assimilated with household waste and follow the household waste circuit.

Objects and materials soiled by blood and/or body fluids must be processed as "Infectious Medical Waste" (IMW).

12. LIMITS OF USE OF INSTRUMENTS

The limits of use of instruments are shown in figures 1 to 13 of this instruction booklet:

- figure 1: residue,
- figure 2: flaking,
- figure 3: cracking,
- figure 4: wearing of laser marking,
- figure 5: blunt, worn cutting parts,
- figure 6: corrosion,
- figure 7: scratches,
- figure 8: defect at the attachment,
- figure 9: metal particles,
- figure 10: break,
- figure 11: dulling,
- figure 12: deformation,
- figure 13: seizing.

Key references:

Complies

Instrument to monitor

At end of life

- figures 16 and 17: reamer holder TFT-AO,
- figure 18: liner head impactor Clicker,
- figure 19: long expansive straight.

14. DATE INSTRUMENTS FIRST PUT ON THE MARKET:

- class I hip instruments: 2022 (MDR),
- class I extremity instruments: 2021 (MDR),
- class Ir instruments: 2023 (MDR),
- class IIa instruments: 1992 (MDD),
- class I hip baskets: 2024 (MDR).

We would like to thank you for your confidence and remain at your disposal for any further information.

13. INSTRUCTIONS FOR ASSEMBLY AND DISASSEMBLY OF INSTRUMENTS

The instructions for assembling and disassembling instruments, when applicable, are presented in figures 14 to 19 of this instruction booklet:

- figure 14: head impactor ICJB E,
- figure 15: depth gauge MV 700,

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